



Power & long life on the go in a battery operated suction unit!

It Performs ...

No need to worry about power when you need it. This high performance unit delivers a vacuum range of 2-21 inHg. The standard easy-to-read gauge allows the care giver to adjust the vacuum.

It Lasts...

The high performance battery allows 1.5 hours of run time, twice that of comparable units. A "Low Battery Light", standard in this unit, informs the caregiver that the unit needs to be re-charged thus preventing interruptions in use.

It's Convenient ...

No need to buy a special canister. The Easy Go Vac uses any standard 800 cc canister and is visible in the Carry Bag while unit is in operation. It is re-chargeable in its convenient, specially designed Carry Bag, which also stores the charger, Suction Tubing and optional Auto DC Power Cord.

It's Durable ...

The ABS plastic cabinet is durable, attractive, easy-to-clean and infection control friendly. There are no open areas on the unit to harbor bacteria and dirt.

It's Protected ...

A dual pump protection system features two safeguards. A float shut-off in the canister to protect overflow and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump.

It's Dependable ...

The Easy Go Vac carries a full 2 year warranty on the unit and a 6 month warranty on the battery.



SPECIFICATIONS

Dimensions (in Carry Bag):

Length 12" (30.3 cm), Width $9^{1/2}$ " (24.1 cm), Height 10³/₁₆" (25.9 cm)

Weight (with Carry Bag and all Accessories): 9.6 lbs. (4.4 kg)

Shipping Weight:

11 lbs. (5 kg)

Electrical:

Aspirator unit: 12 Volts DC Battery Charger/Power Supply: 100-240 Volts AC

Operating Temperature Range:

50°F to 104°F (10°C to 40°Č)

Storage Conditions:

Temperature: $-4^{\circ}F$ to $104^{\circ}F$ ($-20^{\circ}C$ to $40^{\circ}C$); Humidity: Max 95% non-condensing

Vacuum Range:

50 mmHg to 535 mmHg (2 inHg to 21 inHg)

Flow:

14 lpm (unrestricted flow, without suction canister, running on AC power)

Power Consumption:

20 watts

Vacuum Buildup:

0 to 300 mmHg in 3.1 sec. (with 800cc Suction Canister and AC Power Cord connected) 0 to 300 mmHg in 3.5 sec. (with 800cc Suction Canister connected and running on Battery)

Suction Canister:

800cc Disposable with Float Shut-off

ORDERING INFORMATION

Part#	Description
PM65HG	Easy Go Vac Aspirator (includes all items in Ordering Information except Auto DC Power Cord)
1955	Patient Suction Tubing (6ft.)
502519-10	Canister, 800cc (10 pack)
502956	Carry Bag

Part#	Description
PM65HG	Easy Go Vac Aspirator (includes all items in Ordering Information except Auto DC Power Cord)
1955	Patient Suction Tubing (6ft.)
502519-10	Canister, 800cc (10 pack)
502956	Carry Bag



ISO 9001 Certified

Part#	Description
502954	Auto DC Power Cord (optional)
503106	AC Hospital Grade Power Cord
503293	Filter Assembly (Filter, Tubing and Elbow)
503974	Battery Charger/Power Supply

Precision Medical, Inc.

300 Held Drive • Northampton, PA 18067

Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285 Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240 Visit us on the Web: www.precisionmedical.com

503334rev4 6/05 7.5M Printed in USA