



# Disposable Monopolar EMG Needle Electrode - pencil point



At Technomed we understand the importance of patient comfort and signal quality in EMG.

Our new pencil point Monopolar Needle Electrode and cable solutions are designed to give you optimal balance between needle sharpness and control of movement inside the muscle.

The new hub and cable provide a firm grip and optimal precision in needle positioning.

Optional noise reduction cable available.

#### **User Benefits**

- · Clear and reliable signal
- Optimal patient comfort with ease of penetration and repositioning
- · Tip design reduces bleeding
- Precise handling with new ergonomic hub and cable design
- · Approved for stimulation
- Optional 50/60 Hz interference elimination with shielded cable

## Pre-wired

Needle Length	Diameter	Color	Part Number
25 mm	0.35 mm (28 Gauge)	Red	TE/AP2535-335
25 mm	0.45 mm (26 Gauge)	Yellow	TE/AP2545-335
37 mm	0.35 mm (28 Gauge)	Orange	TE/AP3735-335
37 mm	0.45 mm (26 Gauge)	Green	TE/AP3745-335
50 mm	0.45 mm (26 Gauge)	Blue	TE/AP5045-335
75 mm	0.45 mm (26 Gauge)	Violet	TE/AP7545-335

## Detachable

Needle Length	Diameter	Hub Color	Part Number
25 mm	0.35 mm (28 Gauge)	Red	TE/AD2535-008
25 mm	0.45 mm (26 Gauge)	Yellow	TE/AD2545-008
37 mm	0.35 mm (28 Gauge)	Orange	TE/AD3735-008
37 mm	0.45 mm (26 Gauge)	Green	TE/AD3745-008
50 mm	0.45 mm (26 Gauge)	Blue	TE/AD5045-008
75 mm	0.45 mm (26 Gauge)	Violet	TE/AD7545-008

Shelf life for unopened pouches is 36 months 25 individually packed electrodes per box EtO sterilized

Our Reusable Cables for our Detachable Monopolar Needle Electrodes have a perfect fit to our needles and transmit a reliable and qualitative signal.

Cable Type	Wire Length	Part Number
Shielded	1 m	TE/D54-451
Unshielded	0.6 m	TE/D34-231

## Specifications

#### **Properties**

Lead-free	$\checkmark$
Latex-free	$\checkmark$
BPA-free lead wire	
DEHP-free lead wire	$\sqrt{}$
RoHS compliant	$\checkmark$
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802

Pencil point







Your Technomed distributor:





**Technomed Europe** 

Amerikalaan 71 6199 AE Maastricht-Airport The Netherlands T +31 43 408 68 68 info@technomed.nl www.technomed.nl