

Taking you further. Step by step.

# Disposable Hypodermic EMG Needle Electrode

At Technomed we understand the importance of a crystal clear signal, ease of use and patient comfort.

Our Hypodermic Needle is a high quality product developed for EMG guided injections of medication.

Our Hypodermic Needle can be used for multiple insertions without any degradation of its characteristics.

The hub is designed to give a good, secure feel in the hand.

#### **User Benefits**


- Optimal patient comfort with ultra sharp tip and special coating for very smooth penetration and repositioning
- Leak-free and minimum residual volume of medication
- Perfect fit to all slip-tip and Luer lock syringe connections
- Efficiency of use with colour coding, easy to open pouch for quick and hassle free preparation of the procedure

### Available variants

Needle Length	Diameter	Hub Color	Part Number
25 mm	0.30 mm (30 Gauge)	Red	TE/L2530-335
25 mm	0.40 mm (27 Gauge)	Yellow	TE/L2540-335
37 mm	0.40 mm (27 Gauge)	Orange	TE/L3740-335
37 mm	0.45 mm (26 Gauge)	Green	TE/L3745-335
50 mm	0.50 mm (25 Gauge)	Blue	TE/L5050-335
75 mm	0.70 mm (22 Gauge)	Violet	TE/L7570-335

Each Hypodermic Needle has a color coded leadwire (0.75 m) with a DIN 42802 connector

Shelf life for unopened pouches is 36 months 10 individually packed electrodes per box EtO sterilized



## Specifications

### **Properties**

Lead-free electrode	$\checkmark$	
Latex-free electrode	$\checkmark$	
BPA-free lead wire	$\checkmark$	
DEHP-free lead wire	$\checkmark$	
RoHS compliant	$\checkmark$	
Lead wire core	Tin Plated Copper	
Lead wire insulation	PVC	
Connector	DIN 42802	



#### Distributed by:

KEGO USA 1191 Brock McVey Dr. Unit D Lexington, KY 40509

Tel: 800.600.1390 www.kegodiagnostics.com

KEGO Canada 411-419 Consortium Court London, Ontario N6E 2S8

Tel: 866.862.7328 www.kegocorp.com www.kegodiagnostics.ca



D005999 01 Rev C

Taking you further. Step by step.