

Disposable EEG Cup Electrode Silver/Silver Chloride



At Technomed we understand the importance of efficiency, clear signals, and ease of use.

Our disposable Ag/AgCl EEG Cup Electrode is designed to optimize EEG procedures and to enable you to focus on the patient and the EEG recording.

The electrode is designed with flexible wires and its connector is optimized for easy head box insertion and extraction.

User Benefits

- No more cleaning and preparation of EEG electrodes
- Clear and reliable signal with optimized layer of frequency stable Ag/AgCl
- Optimized connector for easy headbox insertion and extraction
- Always new and reliable electrodes at hand
- No risk of cross contamination

Available variants

Wire Length	Part Number
-------------	-------------

1 m	TE/C52-438
1.5 m	TE/C52-638
2 m	TE/C52-838
2.5 m	TE/C52-938

All variants are packaged in sets of 25, 5 x 5 colors: red, green, white, yellow, blue

Specifications

Dimensions

Diameter of cup	10 mm
Height of cup	2.6 mm

Properties

Lead-free	✓
Latex-free	✓
BPA-free lead wire	✓
DEHP-free lead wire	✓
RoHS compliant	✓
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802



USA
Rx ONLY



Your Technomed distributor: