

Reusable EEG Cup Electrode Gold and Silver



At Technomed we understand the importance of a crystal clear signal, ease of use and electrode durability.

Our new Reusable Au and Ag EEG Cup Electrodes are designed to give you reliable signal, and ease of use time after time.

Our new electrode head is designed for easy adjustment to patient contours to avoid electrode lift, and the connector is optimized for easy insertion and extraction.

User Benefits

- Clear and reliable signal with optimized layer of Gold or Silver plating
- Elimination of electrode lift with adjustable electrode neck
- Optimal patient comfort with low profile and soft and flexible anti tangle lead wires
- Easy filling of cups with "curbed" electrode hole
- Optimized connector for easy headbox insertion and extraction
- Long electrode life with Kevlar reinforced wires

Available variants

| Materials | Wire Length | Part Number |
|-------------|-------------|-------------|
| Silver (Ag) | 1 m | TE/C12-434 |
| Silver (Ag) | 1.5 m | TE/C12-634 |
| Silver (Ag) | 2 m | TE/C12-834 |
| Silver (Ag) | 2.5 m | TE/C12-934 |
| Gold (Au) | 1 m | TE/C32-434 |
| Gold (Au) | 1.5 m | TE/C32-634 |
| Gold (Au) | 2 m | TE/C32-834 |
| Gold (Au) | 2.5 m | TE/C32-934 |

All variants are packaged in sets of 12, 2 x 6 colors: red, green, white, yellow, blue, brown

Specifications

Dimensions

| | |
|-----------------|--------|
| Diameter of cup | 10 mm |
| Height of cup | 3.1 mm |

Properties

| | |
|-------------------------|---------------------------------|
| Lead-free | ✓ |
| Latex-free | ✓ |
| BPA-free lead wire | ✓ |
| DEHP-free lead wire | ✓ |
| RoHS compliant | ✓ |
| Lead wire design | Tinsel wire |
| Lead wire core | Kevlar and Silver plated Copper |
| Lead wire insulation | TPU |
| Spring loaded connector | Gold plated Brass/ PVC |
| Connector | DIN 42802 |

Cleaning and disinfection

Cleaning by washing with non-scratching tools and lukewarm water. Disinfecting with alcohol.



| |
|---------|
| USA |
| Rx ONLY |



Your Technomed distributor:


Technomed
 medical accessories

Taking you further. Step by step.

Technomed Europe
 Amerikalaan 71
 6199 AE Maastricht-Airport
 The Netherlands
 T +31 43 408 68 68
 info@technomed.nl
 www.technomed.nl