

Features

- □ Nebulizes 3cc within 7 minutes or less at 7_{LPM} in horizontal or vertical position
- ☐ Larger surface area provided by unique convex cone design
- ☐ Removable "cone" design color coded
- ☐ Scalloped, bottom cuff design
- ☐ One-eighth turn screw on cap


Benefits

- ☐ Allows patient to be in a comfortable position for treatment
- ☐ Maximum capillary action eliminates medication waste
- ☐ Allows for easy, thorough jet cleaning and reassembly
- ☐ Allows the unit to stand alone; simplifies tube attachment
- ☐ Easy to assemble and disassemble, provides a leak- proof seal



Salter 8900 Series Nebulizer* Disposable—single patient use

Built-in quality and reliability... Minimizes patient operational assistance.



Flow Rate	Particle Size (microns)										
(L/min.)	0.14	0.23	0.40	0.62	1.30	2.30	3.20	4.60	6.50	9.20	MMAD
4	3%	5%	10%	21%	25%	26%	6%	2%	1%	2%	1.30
6	4%	9%	15%	26%	25%	16%	2%	1%	2%	1%	1.10
7	4%	11%	16%	27%	22%	14%	2%	1%	1%	1%	1.10
8	4%	12%	17%	30%	20%	12%	2%	0%	1%	1%	1.00
Test Specifications: were performed at 50psi source pressure with 0.9%NaCl. Particle size distribution measured by QCM Cascade impactor.											

Ordering	Inform	nation
OTHER	III OIII	IVVIVIOIT

8900 Series Nebulizer Kits*	Units	Units Part			D (
Disposable—single patient use	(per case) No.	Units (per case)	Part No.	
Full Kit: Nebulizer, anti-drool "T", mouthpiece, 6" reservoir tube, 7' supply tube	50	8900	Aerosol "T" 22mm O.D. x 22mm O.D. x 15mm I.D.	. 50	8917
Same as above with 7' supply tube with female thread-grip connector	50	8900TG	Nebulizer, adult over-the-ear style aerosol mask	50	8903
In-Line Kit: Nebulizer, 22mm O.D. x 22mm O.D. x 18mm I.D."T", 6" reservoir tube, 7' supply tube	50	8915	Nebulizer, adult elastic strap style aerosol mask	50	8923
Same as above with 7' supply tube with female thread-grip connector	50	8915TG	Nebulizer, adult over-the-ear style aerosol mask, 7'supply tube	50	8904
In-Line "T" 22mm O.D. x 22mm O.D. x 18mm I.D.	50	8916	Same as above with female thread grip connector	50	8904TG
Nebulizer, anti-drool "T", mouthpiece, 7' supply tube	50	8911	Nebulizer, adult elastic strap style aerosol mask, 7' supply tube	50	8924
Nebulizer, anti-drool "T", mouthpiece, 6" reservoir tube	50	8912	Nebulizer,pediatric elastic strap style aerosol mask	50	8905
Nebulizer, anti-drool "T", mouthpiece	50	8913	Nebulizer, pediatric elastic strap style aerosol mask, 7' supply tube	50	8906
Nebulizer, anti-drool "T", 6" reservoir tube	50	8914	Nebulizer, pediatric aerosol mask VADS*, with elastic head strap,	50	8907
Nebulizer, 7' supply tube	50	8902	and 7' supply tube		
Nebulizer only	50	8901		ALTER	THR2 ®

8910

1003

50

Anti-Drool "T" 22mm O.D. x 22mm I.D. x 18mm I.D.

Mouthpiece

Printed in U.S.A. *U.S. Patent No. 4,746,067, RE33,717 © Copyright, 1998 Salter Labs® Revised Sept, 2003 SLC-21

Phone: 661-854-3166 1-800-421-0024 (CA only) 1-800-235-4203 (U.S. & Canada) Fax: 661-854-3850 Toll Free Fax: 1-800-628-4690 (U.S. & Canada) International Fax: 661-854-3850

100 W. Sycamore Road - Arvin, CA 93203

International Fax: 661-854-3850 E-mail: salterlabs@us.salterlabs.com www.salterlabs.com

€ 0482

^{*}VADS (Valved Aerosol Delivery System SLC-43)